

CURLIN[®] 6000 Ambulatory Infusion Pump

BY MOOG[®]

Worksheet for Non-Library Programming

TPN THERAPY

Refer to the CURLIN 6000 User Manual for complete programming instructions.

PATIENT:	PUMP SERIAL #:
----------	----------------

PROGRAMMING / SETUP SEQUENCE

NOTE: Use the **↑ (UP)** or **↓ (DOWN) ARROW** keys to move highlighter and the numeric keys for data entry. Use the **YES** and **NO** keys to accept/confirm or reject/change a field entry.

1. Turn the pump on and wait for the self-test to complete, note the power source with power graph on the screen.

2. Highlight **PROGRAM** from the SELECT menu and press **YES**

NOTE: If the following screen displays "NEW PROGRAM", the PREVIOUS RX has not been cleared. Refer to the last section of the "CURLIN 6000 TRAINING & QUICK REFERENCE GUIDE FOR BETWEEN PATIENT USE" to clear all Previous RX prior to new patient use.

3. Highlight **TPN** from the therapy selection menu and press **YES**.

NOTE: If any on-screen alerts appear, the Previous RX has not been cleared. Refer to the last section of the "CURLIN 6000 TRAINING & QUICK REFERENCE GUIDE FOR BETWEEN PATIENT USE" to clear all Previous RX prior to new patient use.

4. With **Unidentified** highlighted, press **YES** (this indicates Manual Programming Mode)

5. **DELAY:** Press **YES** to leave the delayed start feature off, press **NO** to set a delay then press **YES** to ON:

DELAY:	<input type="text"/>
TIME:	<input type="text"/>
DATE:	<input type="text"/>

TIME: Enter time infusion to begin, press **YES**

DATE: Enter date Infusion to begin, press **YES**

ACCEPT?: Press **YES**

6. At **NEXT?** Press **YES** to continue

7. **BAG VOL:** Enter bag volume, press **YES**. (Bag volume should include overflow volume if pump priming and/or a KVO is programmed.)

BAG VOL:	<input type="text"/>
----------	----------------------

NOTE: To force an auto-calculation, enter any 3 of the 4 entries: Vol TBI, INF RATE, UP/DN Ramp or TOT TIME, the unknown field will auto-calculate. Bypass the unknown field by using the **DOWN ARROW** key (typically the INF RATE). Enter the known fields pressing **YES** to confirm each. The unknown field will auto-calculate. Press **YES** to confirm this field.

8. **Vol TBI:** Enter the amount to be infused, press **YES**

Vol TBI:	<input type="text"/>
----------	----------------------



tel: 801.264.1001
toll free: 800.970.2337
www.curlinpump.com

MOOG[®] Enhancing Healthcare.
Enriching Lives.

9. INF RATE: Enter the infusion rate, press YES	Inf RATE: <input type="text"/>
10. UP Ramp: Enter the up ramp time (if any), press YES	UP Ramp: <input type="text"/>
11. DN Ramp: Enter the down ramp time (if any), press YES	DN Ramp: <input type="text"/>
12. TOT TIME: Enter the total time for the infusion, press YES . Continue pressing YES if needed until the KVO field is highlighted.	TOT TIME: <input type="text"/>
13. KVO RATE: Enter KVO (Keep Vein Open) rate, then press YES (0 to 10 ml/hr)	KVO RATE: <input type="text"/>

14. **DONE?:** Press **YES** to accept

When at the **Run to Start Screen**, press the **HELP/OPTIONS** key to set the Lock Level (security level), and to change the Air-In-Line sensitivity and the Down Occlusion sensitivity from the default settings when the therapy requires these changes:

Lock: Press **YES** to accept or **NO** to change Lock setting, press **YES** to accept new setting.

Lock	<input type="text"/>
------	----------------------

DN Occlu: Press **YES** to accept (default setting is LOW) or press **NO** to change then press **YES** to accept.

DN Occlu:	<input type="text"/>
-----------	----------------------

AIR SENS: Press **YES** to accept (default setting is 0.5 mL) or **NO** to change then press **YES** to accept.

AIR SENS:	<input type="text"/>
-----------	----------------------

Press **YES** at each field to advance to the next field. Press **YES** to **ACCEPT OPT** at the end of the menu.

PROGRAMMED BY:	DATE:
CHECKED BY:	DATE:



tel: 801.264.1001
 toll free: 800.970.2337
 www.curlinpump.com

MOOG® Enhancing Healthcare.
 Enriching Lives.

CURLIN® is a registered trademark of Curlin Medical, Inc. Moog® is a registered trademark of Moog, Inc. © 2018 Moog, Inc. Moog Medical, 4314 Zevex Park Lane, Salt Lake City, UT 84123 USA. LS-61476 Rev. A